

## SE-1200 Express SE-1200 Express Basic 12-Channel ECG



The Edan SE-1200 ECG machines boast exceptional cost performance with both accurate clinical performance as well as user-friendly design. These machines feature Glasgow Interpretation making your workflow more accurate and efficient. The SE-1200 ECG machines can store up to 800 ECGs.

### What's The Difference?

SE-1200 Express Basic

Short-cut keys only

SE-1200 Express

Touch screen and short cut keys



Wireless AP



One-Touch Operation

# User-Friendly Design

## Highlighted Features

### SE-1515 PC Software



- Compare previous ECGs
- Increase confidence when diagnosing
- Easily export ECG reports

### Data Management

On board PDF creator to transmit data to USB

Optional SE-1515 PC Software

\*\*must be purchased or utilized to send wireless reports\*\*

### Exceptional Clinical Performance

Weak signal detection with 24 bits A/D converter technology

Anti-noise technology with comprehensive filters

High CMRR over 115dB

120s Real-time ECG waveforms freezing and review

Automatic measurements and interpretation tested by CSE/AHA/MIT database

Extended printing when arrhythmia is detected

Integrated advanced ECG analysis: Vector ECG



Built-in thermal recorder



8-inch LCD screen



USB support



Alphanumeric keyboard



Barcode scanner support



One-touch operation



Shortcut keys



PDF, SCP, FDA-XML, DICOM

DATA transmitting to PC via LAN/WIFI

PDF, SCP, FDA-XML and DICOM data export options