



Better Cardiac Care



PCA-500

PERSONAL ECG

PERSONAL 12-LEAD ECG


Now you can get hospital-quality ECG anywhere, anytime, in the comfort of your home.

PCA 500


The streamlined ECG workflow is time-efficient, cost-effective and optimized for better quality of care.




Conventional 12-lead ECG workflow



Workflow of PCA by QT Medical




CALL FOR MORE INFO
888.602.PROS




Hospital-quality 12-lead ECG with
improved Efficiency, Accuracy, Consistency –

NO need for a technician



Conventional



PCA 500



Efficiency

Performing an ECG with the PCA 500 takes only 3 minutes.
Comparing to conventional ECG, using the PCA 500 Electrode
can save up to 70% time.



Accuracy

The intuitive design of the system allows patients with no
training to perform a full 12-lead ECG on themselves without a
technician.



Consistency

The PCA 500 Electrode enables consistent leads placement,
thus subtle changes in ECG over time can be detected.


 **MDPRO**
PCA-500 PERSONAL ECG

CALL FOR MORE INFO
888.602.PROS



Complete ECG solution Anywhere, Anytime

PCA 500 is wireless,
digital and high quality. It
utilizes smart mobile devices
and manages in the cloud.
There are many places with
needs for accurate ECG but
could not get it before— now
PCA 500 is the perfect
solution!



Telemedicine



Ambulances



ER's &
Urgent Care



CRO/ Clinical
Trials



Employee
Health



Athlete
Screening



Remote Area



Post-Op,
CABG, PCI



Nursing
Homes, SNF



Jails, Prisons



Airlines




Home Care

 **MDPRO**
PCA-500 PERSONAL ECG

CALL FOR MORE INFO
888.602.PROS




PCA 500 - Personal Cardiac Assistant



Compact and portable resting 12-lead ECG system

- Build-in 700 mAh rechargeable Li-ion battery
- 16 hours of continuous use
- Wireless connectivity via Bluetooth
- Leads-off indicator
- Dimension: 7.2 x 6.8 x 1.8 cm
- Weight: 67 g



Patented pre-positioned electrodes

- Simplifies lead placement procedure
- 70% faster setup time
- AHA & IEC compliant
- Proven usable by patients with no training



QT ECG Apps & cloud management system

- Apps available for iOS and Android phones & tablets
- 1,000Hz sampling rate suitable for pediatrics and adults
- Adjustable gain and speed function
- Serial comparison of different ECG's by superimosing the waveforms
- All documents and encrypted data are transferred through a secure connection to the cloud compliant with HIPAA & GDPR



PCA-500 PERSONAL ECG

CALL FOR MORE INFO
888.602.PROS